



PATENT

Docket No.: I9000.0058/P058

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re Patent Application of:
Marit Nilsen-Hamilton

BEST AVAILABLE COPY

Serial No.: 10/809,886

Group Art Unit: 1635

Filed: March 26, 2004

Examiner: Kimberly Chong

For: ALLOSTERIC PROBES AND
METHODS

Assistant Commissioner for Patents
Washington, D.C. 20231

DECLARATION OF MARIT NILSEN-HAMILTON UNDER
37 CFR 1.131

Dear Sir:

I, Marit Nilsen-Hamilton, declare and state as follows:

1. I am of legal age, and under no disability that prevents me from attesting to the following statements and information, which are based on my personal knowledge and observations or on my best information and belief.

2. I reside at 1111 N. Hyland Avenue, Ames, Iowa 50014.

3. I am the sole inventor of the above-identified U.S. patent application (the "'886 application"), filed on March 26, 2004, as evidenced by the attached executed Declaration (Exhibit A).



4. I reviewed and understand the above-identified U.S. patent application including the currently pending claims and amendments (the "Claimed Invention").

5. I conceived of and reduced to practice the invention covered by the Claimed Invention prior to April 24, 2002. My reduction to practice of the Claimed Invention is evidenced by Exhibits B-D, which are summaries documenting work I directed and supervised. Specifically, Exhibits B-D are photocopies of handwritten notes, and electronic summaries of the handwritten notes, both created by my student who conducted experiments under my direction and supervision. Each of the experiments summarized in Exhibits B-D was conducted by my student under my direction and supervision prior to April 24, 2002, although the electronic summaries were printed after April 24, 2002. The actual dates that the experiments were conducted have been blanked out, as has any description not relevant to the conception of the Claimed Invention; however, the dates are all prior to April 24, 2002.

6. Exhibits B-D shows the results of an experiments involving CLAMP 609, which is a probe for binding a plurality of targets comprising an allosteric regulator linked to at least one regulated aptamer. Specifically, as noted in Exhibit D, CLAMP 609 includes an allosteric regulator, which is a nucleic acid molecule, linked to a regulated aptamer in cis configuration.

Serial No.: 10/809,886

Docket No.: 19000.0058/P058

7. The experimental results described in Exhibits B-D show that binding the allosteric regulator with a first target enhances the binding of the at least one regulated aptamer to at least one second target. Specifically, as shown in Exhibits B-D, when neomycin binds the allosteric regulator of the CLAMP 609, the regulated aptamer is increasingly able to bind ATP.

8. Thus, Exhibits B-D show that, prior to April 24, 2002, the Claimed Invention was reduced to practice under my direction and supervision.


All statements made herein of my own knowledge are true, and all statements made on information and belief are believed to be true. All statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the above-identified patent.

Date: 9/27/06By: Marit Nilsen-Hamilton
Marit Nilsen-Hamilton

Serial No.: 10/809,886

Docket No.: 19000.0058/P058

EXHIBIT A



**Declaration for Patent Application
English Language Declaration**

Attorney Docket No. | 19000.0058/P058
First Named Inventor | Marit Nilsen-Hamilton

COMPLETE IF KNOWN:

☐ Submitted with initial filing
☒ Submitted after initial filing (surcharge required 37 CFR 1.16(e))

Application No. | 10/809,886-Conf. #7776
Filing Date | March 26, 2004
Art Unit | 1645
Examiner | Not Yet Assigned

As a below named inventor, I hereby declare that:

My residence, mailing address and citizenship are as stated below next to my name.

I believe I am the original, first and sole inventor (if only one name is listed below) or an original, first and joint inventor (if plural names are listed below) of the subject matter which is claimed and for which a patent is sought on the invention entitled:

ALLOSTERIC PROBES AND METHODS

the specification of which

☐ is attached hereto
OR

☒ was filed on 03/26/2004

as United States Application No. or PCT International Application No. 10/809,886
and was amended on (if applicable).

I hereby state that I have reviewed and understand the contents of the above-identified specification, including the claims, as amended by any amendment referred to above.

I acknowledge the duty to disclose information which is material to patentability as defined in 37 CFR 1.56, including for continuation-in-part applications, material information which became available between the filing date of the prior application and the National or PCT International filing date of the continuation-in-part application.

I hereby claim foreign priority benefits under 35 U.S.C. 119(a)-(d) or (f), or 365(b) of any foreign applications(s) for patent, inventor's or plant breeder's rights certificate(s), or 365(a) of any PCT international application which designated at least one country other than the United States of America, listed below and have also identified below, by checking the box, any foreign application for patent, inventor's or plant breeder's right certificate(s), or any PCT international application having a filing date before that of the application on which priority is claimed.

Prior Foreign Application(s)

			Priority Not Claimed	Certified Copy Attached	
				YES	NO
(Number)	(Country)	(Filing Date)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
(Number)	(Country)	(Filing Date)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
(Number)	(Country)	(Filing Date)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

☐ Additional prior foreign applications are listed on a supplemental data sheet attached hereto.

I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under 18 U.S.C. 1001 and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

Full name of sole or first inventor Marit Nilsen-Hamilton	
Sole or first inventor's signature <i>Marit Nilsen Hamilton</i>	Date 9/2/04
Residence Ames, Iowa	
Citizenship Norway	
Mailing Address 1111 N. Hyland Ave. Ames, Iowa 50014-4007	

STATE OF IOWA)
) ss:
COUNTY OF STORY.)

On this 2nd day of September, 2004, before me personally came Marit Nilsen-Hamilton, to me known to be the individual described in and who executed the foregoing instrument, and acknowledged execution of the same.

Serial No.: 10/809,886

Docket No.: 19000.0058/P058

EXHIBIT B

TITLE: Column Affinity - CLAMP 609 on ATP Agarose in Neomycin Solution

DATE: _____

PURPOSE: To test the effect of free neomycin on Neo/ATP-CLAMP 609 binding

EXP.# BK036

REAGENTS

DEPC treated water

Neo/ATP-CLAMP 609 (~ 0.5pmol/ul) BK035 (gcuaauuacgacucacuaauaggccugggcgagaaguuuaggccuu
ggguugggaagaaacuguggcacuucggugccagcaaccc)

Neomycin Sepharose from BK001

ATP Agarose from Sigma No. A2767

MicroBioSpin Column from BioRad No. 732-6204

Neomycin from Sigma No. N-1876 Lot 90K0854 MW 908.9

Adenosine 5'Triphosphate (ATP) from Sigma No. A-3377 MW 551.1

Binding Buffer: 50mM Tris, 250mM NaCl, 5mM MgCl₂, pH 7.6

Binding Buffer (ATP columns): 1mM neomycin in Binding Buffer above

Elution Buffers: 10mM Neomycin or 10 mM ATP in Binding Buffer

Fisher Scientific ScintiVerse Scintanalyzed Cocktail No. SX18-4

Packard TR1600 Liquid Scintillation Analyzer (1238 Molec Biol) - Protocol 3

STEP	SOLUTION	VOL (ml)	TEMP (°C)	Time Period (l)
Make Column	Add 100ul either neomycin sepharose or ATP agarose to each column.	100ul	23	5min
Wash Column	Equilibrate column by washing 10x with 100ul Binding Buffer, collecting 1ml.	1x1ml/fx	23	5min
Prepare sample	Add ~ 1pmol ³² P-CLAMP609 to DEPC water for each sample and preincubate.	20ul	85 23	1min 1min
Apply Sample	Apply sample to column.	18ul	23	30sec
Wash	Wash column 20x with 100ul Binding Buffer, collecting each 1ml.	2x1ml/fx	23	10min
Elution	Elute 10x with 100ul Elution Buffer, collecting each 100ul.	10x100ul/fx	23	10min
**Column Wash	Wash column 10X with 100ul Elution Buffer, then 10x with 100ul Binding Buffer, collecting each 1ml.	2x1ml/fx	23	10min
Analyze Fractions	Add 20ul each fraction to 5ml ScintiVerse and assay on liquid scintillator. Reserve 2ul of sample for total sample cpm count.	20ul in 5ml	23	112min

**Column Washes only for ATP agarose column

DISCUSSION:

The neomycin in the binding buffer apparently increased the effective ATP binding of CLAMP 609. In the first trial, 90% of the TCA precipitable counts bound to the neomycin sepharose and 17% were eluted from the ATP agarose (with only about 50% of the total radioactivity recovered from the column).** In the second trial, again about 90% bound to the neomycin sepharose and 35% bound to the ATP agarose, with much less radioactivity remaining in the column (about 80% recovered). This is a great improvement over the initial CLAMP 609 binding assays (BK033 & BK034) where less than 10% was shown to bind to the ATP agarose. Next this difference in binding must be directly compared in a single experiment with and without neomycin in the binding buffer.

**Note, the elution in the first trial with ATP was distributed over the entire elution series and was probat background. However, in the second trial there was a definite peak of CLAMP binding

INIT____DATE:____PRINTED: May 13, 200:

Expt BK038
printed 13-May-02
CLAMP 609

Neomycin Sepharose		dilution			Acid precipitable		
sample	value	value - bkg	factor	total cpm	%input	cpm	% total
input	14712.5	14672	9	132048		129407.0	100%
wash	387.5	347	50	17350	13.14%		0%
wash	40.5	0	50	0	0.00%		
1	9017.5	8977	5	44885	33.99%	44885.0	34.89%
2	7018.5	6978	5	34890	28.42%	34890.0	28.96%
3	3138	3095.5	5	15477.5	11.72%	15477.5	11.96%
4	1895	1854.5	5	9272.5	7.02%	9272.5	7.17%
5	1284.5	1244	5	6220	4.71%	6220.0	4.81%
6	874.5	834	5	4170	3.16%	4170.0	3.22%
7	573.5	533	5	2665	2.02%	2665.0	2.08%
8	413.5	373	5	1865	1.41%	1865.0	1.44%
9	309	288.5	5	1342.5	1.02%	1342.5	1.04%
10	298.5	256	5	1280	0.97%	1280.0	0.99%

BKG 40.5 Eluted
% acid Insoluble 106% 122068 94.33%
98.0%

ATP Agarose		dilution			Acid precipitable		
sample	value	value - bkg	factor	total cpm	%input	cpm	% total
input	21260.5	21220	9	190980		187160.4	100%
wash	1026.5	986	50	49300	25.81%		0%
wash	173.5	133	50	8650	3.48%		
1	289	248.5	5	1242.5	0.85%	1242.5	0.68%
2	599.5	559	5	2795	1.46%	2795	1.49%
3	853.5	813	5	4065	2.13%	4065	2.17%
4	780.5	740	5	3700	1.94%	3700	1.98%
5	867	826.5	5	3132.5	1.64%	3132.5	1.67%
6	728.5	688	5	3440	1.80%	3440	1.84%
7	719	678.5	5	3392.5	1.78%	3392.5	1.81%
8	886	845.5	5	3227.5	1.69%	3227.5	1.72%
9	683	642.5	5	3212.5	1.68%	3212.5	1.72%
10	592.5	552	5	2780	1.45%	2760	1.47%

BKG 40.5 Eluted
% acid Insoluble 48% 30968 16.55%
98.0%


Neomycin Sepharose			dilution		Acid precipitable		
sample	value	value - bkg	factor	total cpm	%Input	cpm	% total
Input	8283	8283	9	74547		73056.1	100%
wash	67.5	67.5	50	3375	4.53%		
wash	23.5	23.5	50	1175	1.58%		
1	2544.5	2544.5	5	12722.5	17.07%	12722.5	17.41%
2	2970.5	2970.5	5	14852.5	19.92%	14852.5	20.33%
3	2268	2268	5	11340	15.21%	11340.0	15.52%
4	1453	1453	5	7265	9.75%	7265.0	9.94%
5	1196	1196	5	5980	8.02%	5980.0	8.19%
6	977.5	977.5	5	4887.5	6.56%	4887.5	6.69%
7	608	608	5	3040	4.08%	3040.0	4.16%
8	499	499	5	2495	3.35%	2495.0	3.42%
9	334.5	334.5	5	1672.5	2.24%	1672.5	2.29%
10	262.5	262.5	5	1312.5	1.76%	1312.5	1.80%

BKG 23.5 Eluted
% acid Insoluble 94% 65568 89.75%
98.0%

ATP Agarose		dilution		Acid precipitable			
sample	value	value - bkg	factor	total cpm	%input	cpm	% total
input	8184	8127.5	9	73147.5		71684.55	100%
wash	596	539.5	50	26975	36.88%		0%
wash	171.5	115	50	5750	7.86%		
1	587.5	531	5	2655	3.63%	2655	3.70%
2	2057	2000.5	5	10002.5	13.67%	10002.5	13.95%
3	1423	1366.5	5	6832.5	9.34%	6832.5	9.53%
4	686.5	630	5	3150	4.31%	3150	4.39%
5	350.5	294	5	1470	2.01%	1470	2.05%
6	162.5	106	5	530	0.72%	530	0.74%
7	93.5	37	5	185	0.25%	185	0.26%
8	67.5	11	5	55	0.08%	55	0.08%
9	56.5	0	5	0	0.00%	0	0.00%
10	57	0.5	5	2.5	0.00%	2.5	0.00%

BKG 58.5 Eluted
% acid insoluble 79% 24883 34.71%
98.0%

BK038: CLAMP 609 on affinity gels



—■— Neomycin Sepharose 1 —▲— ATP Agarose 1
- - ■ - - Neomycin Sepharose 2 - - ▲ - - ATP Agarose 2.

20136

[Protocol] #: 3

Name: 32PMNH

12.21

S#

CLAMP 609 binding to ATP gel w/ "unlocking" by neomycin

20136

CLAMP 609 (840³³)

Neomycin Sepharose (new) for control; ATP agarose (used)

Binding Buffer (Neo): 50 mM Tris, 250 mM NaCl, 5 mM MgCl₂, pH 7.6

Binding Buffer (ATP): 1 mM neomycin in above

Elution Buffer (Neo): 10 mM neomycin in

Elution Buffer (ATP): 10 mM ATP in

} old solns → need to make fresh

Wash column 10x w/ 100 µl Binding Buffer.

Add 4 µl ^{PNA} → 20 µl sample. Remove 2 µl → 20 µl for total cpm count.

Incubate remaining sample 1 min @ 85°C, cool 1 min, add to column.

Wash 20x w/ 100 µl Binding Buffer, collecting each 1 ml.

Elute 10x w/ 100 µl Elution Buffer, collecting each 100 µl.

Assay 20 µl / fx in 5 ml Scintivial

* Forced drip-quick flow

* Columns still hot after elution

Samples: 1 Control 20 µl/20 µl 1/4 sample

2-3 Washes 20 µl/100 µl

4-13 Elution 20 µl/100 µl

14 Control 20 µl/20 µl 1/4 sample

15-16 Washes 20 µl/100 µl

17-26 Elution 20 µl/100 µl

27-28 Column wash 20 µl/100 µl

} Neomycin sepharose

} ATP agarose

⇒ still low elution from ATP column

170-5V

Protocol #: 3

Name: 32PMNH

2 12.21

S#

CLAMP 609 Binding to ATP gel w/ "unlocking" by neomycin (repeat)

36(v)

CLAMP 609 (24022)

Neomycin sepharose (neo) for control $\frac{1}{2}$ ATP agarose (used)(Neo) Binding Buffer: 50 mM Tris, 250 mM NaCl, 5 mM MgCl₂, pH 7.6(ATP) ~~Binding Buffer~~: 1 mM neomycin in \rightarrow Elution Buffer (Neo): 10 mM neomycin in \rightarrow Elution Buffer (ATP): 10 mM ATP in \rightarrow Wash column 10x w/ 100 μ l Binding Buffer.Add 2 μ l RNA \rightarrow 20 μ l sample. Remove 2 μ l \rightarrow 20 μ l for total cpm.

Incubate remaining sample 1 min @ 85°C, cool 1 min, add to column.

Wash 20x w/ 100 μ l Binding Buffer, collecting each 1 μ l.Elute 10x w/ 100 μ l Elution Buffer, collecting each 100 μ l.Assay 20 μ l / μ l in 5 μ l Scintiverse.

* Slower drip rate

* Columns only slightly hot after elution (\sim 300 cpm)

Samples: 1	Control cpm	20 μ l / 20 μ l	$\frac{1}{4}$ sample	} Neomycin sepharose
2-3	Washes	20 μ l / 100 μ l		
4-13	Elution	20 μ l / 100 μ l		
14	Control cpm	20 μ l / 20 μ l	$\frac{1}{4}$ sample	} ATP agarose
15-16	Column / rawashes	20 μ l / 100 μ l		
17-18	Washes	20 μ l / 100 μ l		
19-28	Elution	20 μ l / 100 μ l		
29-30	Column Washes	20 μ l / 100 μ l		

No 334

Protocol #: 3 Name: 32PMNH

Region A: LL-UL= 5.0-1700 Ler= 0 Bkg= 0.00 %2 Sigma=0.10
Region B: LL-UL=50.0-1700 Ler= 0 Bkg= 0.00 %2 Sigma=0.10
Region C: LL-UL= 0.0- 0.0 Ler= 0 Bkg= 0.00 %2 Sigma=0.00
Time = 2.00 QIP = SIS

13:21

S#	TIME	CPMA	CPMB	B:2S%	SIS	FLAG
1	2.00	14712.5	11394.5	1.32	627.11	
2	2.00	387.50	304.00	8.11	628.76	
3	2.00	40.50	25.50	28.01	381.93	
4	2.00	9017.50	8926.50	1.70	618.68	
5	2.00	7018.50	5178.00	1.97	553.75	
6	2.00	3136.00	2218.00	3.00	449.81	
7	2.00	1895.00	1422.00	3.75	562.68	
8	2.00	1284.50	872.50	4.79	428.33	
9	2.00	874.50	559.00	5.98	405.60	
10	2.00	573.50	438.50	6.77	583.56	
11	2.00	413.50	312.00	8.01	564.44	
12	2.00	309.00	224.50	9.44	522.36	
13	2.00	296.50	202.50	9.94	459.84	
14	2.00	21260.5	17039.0	1.08	679.43	
15	2.00	1026.50	780.50	5.08	581.78	
16	2.00	173.50	114.50	13.22	444.05	
17	2.00	289.00	198.50	10.09	476.83	
18	2.00	599.50	428.00	8.84	544.94	
19	2.00	853.50	621.50	5.67	561.07	
20	2.00	780.50	594.50	5.80	602.36	
21	2.00	667.00	437.00	8.77	422.81	
22	2.00	728.50	528.00	6.15	483.34	
23	2.00	719.00	521.00	6.20	516.21	
24	2.00	686.00	506.50	8.28	557.88	
25	2.00	683.00	562.00	5.97	690.32	
26	2.00	592.50	494.00	8.38	700.70	
27	2.00	623.00	437.00	8.77	535.24	
28	2.00	108.50	81.00	15.71	495.06	

Protocol #: 3

Name: 32PMNH

Region A: LL-UL= 5.0-1700

Lor=

0

Bkg= 0.00

%2 Sigma=0.10

10:21

Region B: LL-UL=50.0-1700

Lor=

0

Bkg= 0.00

%2 Sigma=0.10

Region C: LL-UL= 0.0- 0.0

Lor=

0

Bkg= 0.00

%2 Sigma=0.10

Time = 2.00

QIP = SIS

S#	TIME	CPMA	CPMB	B:2S%	SIS	FLAG
1	2.00	8283.00	6509.50	1.75	616.07	
2	2.00	67.50	44.00	21.32	570.14	
3	2.00	23.50	10.50	43.64	239.08	
4	2.00	2544.50	1837.00	3.30	482.79	
5	2.00	2970.50	2130.00	3.06	532.72	
6	2.00	2268.00	1735.50	3.39	584.36	
7	2.00	1453.00	1134.50	4.20	640.86	
8	2.00	1196.00	918.50	4.67	635.13	
9	2.00	977.50	791.50	5.03	711.35	
10	2.00	608.00	403.00	7.04	416.83	
11	2.00	499.00	384.00	7.22	543.59	
12	2.00	334.50	228.00	9.37	429.98	
13	2.00	262.50	185.50	10.38	484.77	
14	2.00	8184.00	8188.50	1.80	555.04	
15	2.00	21.50	9.50	45.88	247.09	
16	2.00	22.50	10.50	43.64	465.13	
17	2.00	596.00	442.50	6.72	538.64	
18	2.00	171.50	122.00	12.80	478.34	
19	2.00	587.50	409.00	6.99	501.28	
20	2.00	2057.00	1580.00	3.56	591.29	
21	2.00	1423.00	1124.00	4.22	635.41	
22	2.00	686.50	565.50	5.95	693.13	
23	2.00	350.50	288.50	8.33	694.04	
24	2.00	162.50	120.00	12.91	507.30	
25	2.00	93.50	67.50	17.21	560.65	
26	2.00	67.50	43.00	21.57	490.62	
27	2.00	56.50	41.50	21.95	501.46	
28	2.00	57.00	38.50	22.79	496.97	
29	2.00	77.50	48.50	20.31	452.02	
30	2.00	39.50	23.50	29.17	291.04	

Serial No.: 10/809,886

Docket No.: I9000.0058/P058

EXHIBIT C

TITLE: Column Affinity - Neomycin Unlocking of CLAMP 609 on ATP agarose	DATE: 09/28/2006
PURPOSE: To determine if neomycin binding can increase ATP binding ability of Neo/ATP-CLAMP 609; repeated three times	EXP.# BK038

REAGENTS

DEPC treated water
³²P-Neo-uu-ATP-CLAMP 609 (~ 1.3pmol/ul)/BK037
 (gcuuaauacgacucacuaauaggccugggcgagaaguuuaggcc uu
 gggguugggaagaaacuguggcacuuucggugccagcaaccc)
 ATP Agarose from Sigma No. A2787
 MicroBioSpin Column from BioRad No. 732-6204
 Neomycin from Sigma No. N-1876 Lot 90K0854 MW 908.9
 Adenosine 5'Triphosphate (ATP) from Sigma No. A-3377 MW 551.1
 Binding Buffer: 50mM Tris, 250mM NaCl, 5mM MgCl₂, pH 7.6
 Elution Buffer: 10mM ATP in Binding Buffer
 Fisher Scientific ScintiVerse Scintanalyzer Cocktail No. SX18-4 Lot 994498
 Packard TR1600 Liquid Scintillation Analyzer (1238 Molec Biol) - Protocol 3

STEP	SOLUTION		VOL (ml)	TEMP (°C)	Time Period
Make Column	Layer 100ul of ATP agarose into micro biospin column.		100ul	23	30sec
Wash Column	Equilibrate column by washing 10x with 100ul Binding Buffer, collecting 1ml.	± 1 mM neomycin	1ml	23	5min
Prepare Samples	~ 3pmol ³² P-CLAMP609 Binding buffer	± 1 mM neomycin	20ul	85 23	1min 1min
Apply Sample	Apply sample to column.		18ul	23	30sec
Wash	Wash 20x with 100ul Binding Buffer, collecting each 1ml.	± 1 mM neomycin	2x1ml/fx	23	5min
Elution	Elute 10x with 100ul Elution Buffer, collecting each 1ml.	± 1 mM neomycin	10x100ul/fx	23	10min
Column Washes	Wash column with 1ml Elution Buffer, then with 1ml Binding Buffer, collecting each 1ml.		2x1ml/fx	23	5min
Analyze Fractions	Add 20ul each fraction to 5ml ScintiVerse and assay on liquid scintillator. Reserve 2ul of sample for total sample cpm count.		20ul in 5ml	23	64min

DISCUSSION:

The addition of neomycin to the solutions increased binding to the ATP column from 9% to 28% of TCA precipitable counts in the first trial, from 25% to 59% in the second trial, and from 9% to 25% in the third trial. This suggests that neomycin binding is required in order for the ATP aptamer end of the CLAMP to be fully available for binding. However, the ATP binding is still very low compared to the neomycin binding attained, so other factors probably still reduce this binding, such as incomplete transcription or remaining interference within the CLAMP construct.

CLAMP 609 on ATP Agarose		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	11184	11142.5	9	100282.5				93282.7	100%
wash	1360.5	1339	50	86950	86.8%				0%
wash	189	177.5	50	8875	8.8%				
1	372	350.5	5	1752.5	1.7%	1752.5	1.9%		
2	489.5	468	5	2340	2.3%	2340.0	2.5%		
3	348	327.5	5	1637.5	1.6%	1637.5	1.8%		
4	260.5	239	5	1195	1.2%	1195.0	1.3%		
5	120	98.5	5	492.5	0.5%	492.5	0.5%		
6	72.5	51	5	255	0.3%	255.0	0.3%		
7	53.5	32	5	160	0.2%	160.0	0.2%		
8	37	15.5	5	77.5	0.1%	77.5	0.1%		
9	39	17.5	5	87.5	0.1%	87.5	0.1%		
10	36.5	15	5	75	0.1%	75.0	0.1%		

BKG 21.5 Eluted
% acid insoluble 93.0% 83.68% 8073 9%

CLAMP 609 on ATP Agarose - with neomycin unlocking		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	11832	11811.5	9	106303.5				98862.28	100%
wash	1274.5	1254	50	82700	59%				0%
wash	249	228.5	50	11425	10.75%				
1	555	534.5	5	2672.5	2.51%	2672.5	2.70%		
2	1571.5	1551	5	7755	7.30%	7755	7.84%		
3	1590	1569.5	5	7847.5	7.38%	7847.5	7.94%		
4	1090	1069.5	5	5347.5	5.03%	5347.5	5.41%		
5	462	441.5	5	2207.5	2.08%	2207.5	2.23%		
6	152	131.5	5	657.5	0.62%	657.5	0.67%		
7	131.5	111	5	555	0.52%	555	0.56%		
8	64	43.5	5	217.5	0.20%	217.5	0.22%		
9	55.5	35	5	175	0.16%	175	0.18%		
10	55	34.5	5	172.5	0.16%	172.5	0.17%		

BKG 20.5 Eluted
% acid insoluble 93.0% 95.70% 27608 28%


CLAMP 609 on ATP Agarose		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	5477	5444.5	9	49000.5				45570.47	100%
wash	1013	980.5	50	49025	100%				0%
wash	187	154.5	50	7725	15.77%				
1	520.5	488	5	2440	4.98%	2440	5.35%		
2	829	796.5	5	3982.5	8.13%	3982.5	8.74%		
3	602	569.5	5	2847.5	5.81%	2847.5	6.25%		
4	294.5	262	5	1310	2.67%	1310	2.87%		
5	152.5	120	5	600	1.22%	600	1.32%		
6	86	53.5	5	267.5	0.55%	267.5	0.59%		
7	51.5	19	5	95	0.19%	95	0.21%		
8	40	7.5	5	37.5	0.08%	37.5	0.08%		
9	30.5	-2	5	-10	-0.02%	-10	-0.02%		
10	30	-2.5	5	-12.5	-0.03%	-12.5	-0.03%		

BKG 32.5 Eluted
% acid insoluble 93.0% 139.40% 11558 25%


CLAMP 609 on ATP Agarose - with neomycin unlocking		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	4437.5	4401	9	39809				38836.37	100%
wash	859.5	823	50	41150	104%				0%
wash	197.5	161	50	8050	20.32%				
1	379.5	343	5	1715	4.33%	1715	4.66%		
2	1059	1022.5	5	5112.5	12.91%	5112.5	13.68%		
3	1211	1174.5	5	5872.5	14.83%	5872.5	15.94%		
4	1002.5	968	5	4830	12.19%	4830	13.11%		
5	448.5	410	5	2050	5.18%	2050	5.57%		
6	254.5	218	5	1090	2.75%	1090	2.96%		
7	168.5	130	5	650	1.64%	650	1.76%		
8	113.5	77	5	385	0.97%	385	1.05%		
9	62.5	26	5	130	0.33%	130	0.35%		
10	53	16.5	5	82.5	0.21%	82.5	0.22%		

BKG 36.5 Eluted
% acid insoluble 93.0% 179.55% 21918 59%

BK038: CLAMP 609 w/ Neomycin Unlocking (set 1)



BK038: CLAMP 609 w/ Neomycin Unlocking (set 2)




CLAMP 609 on ATP Agarose		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	8881	8859	9	79704				74124.72	100%
wash	1114	1089	50	54450	68%				0%
wash	218.5	193.5	50	9675	12.14%				
1	303	278	5	1390	1.74%	1390	1.88%		
2	596	571	5	2855	3.58%	2855	3.85%		
3	311	286	5	1430	1.79%	1430	1.93%		
4	142.5	117.5	5	587.5	0.74%	587.5	0.79%		
5	91	66	5	330	0.41%	330	0.45%		
6	40.5	15.5	5	77.5	0.10%	77.5	0.10%		
7	41	18	5	80	0.10%	80	0.11%		
8	32	7	5	35	0.04%	35	0.05%		
9	25.5	0.5	5	2.5	0.00%	2.5	0.00%		
10	42	17	5	85	0.11%	85	0.11%		

BKG 25 Eluted
% acid insoluble 93.0% 89.08% 6873 9%

CLAMP 609 on ATP Agarose - with neomycin unlocking		dilution		total cpm		%input		Acid precipitable	
sample	value	value - bkg	factor					cpm	% total
Input	1554	1507	9	13583				12613.59	100%
wash	114.5	67.5	50	3375	25%				0%
wash	26.5	-20.5	50	-1025	-7.56%				
1	52.5	5.5	5	27.5	0.20%	27.5	0.22%		
2	68.5	19.5	5	97.5	0.72%	97.5	0.77%		
3	98	51	5	255	1.88%	255	2.02%		
4	128	79	5	395	2.91%	395	3.13%		
5	130	83	5	415	3.06%	415	3.29%		
6	132.5	85.5	5	427.5	3.15%	427.5	3.39%		
7	132	85	5	425	3.13%	425	3.37%		
8	126.5	79.5	5	397.5	2.93%	397.5	3.15%		
9	115.5	68.5	5	342.5	2.53%	342.5	2.72%		
10	118	71	5	355	2.62%	355	2.81%		

BKG 47 Eluted
% acid insoluble 93.0% 40.46% 3138 25%

BK038: CLAMP 609 w/ Neomycin Unlocking (set 3)



"Unlocking" effect of neomycin on CLAMP609 ATP binding

CLAMP609 (BKO) 37

ATP agarose (used)

Binding Buffer (control): 50 mM Tris, 250 mM NaCl, 5 mM MgCl₂, pH 7.6

Binding Buffer (unlocking): 1 mM neomycin, 50 mM Tris, 250 mM NaCl, 5 mM MgCl₂, pH 7.6

Election Buffer: 10 mM ATP, 50 mM Tris, 250 mM NaCl, 5 mM MgCl₂, pH 7.6

Election Buffer (unlocking): above plus 1 mM neomycin

Wash column 10x w/ 100 µl Election Buffer, collecting 1 ml.

Wash column 10x w/ 100 µl Binding Buffer, collecting 1 ml.

Add 2 µl PVA → 20 µl sample (B.S. dilutes). Remove 2 µl → 20 µl for total spn.

Incubate remaining sample 1 min @ 85°C, cool 1 min, add to column.

Wash 20x w/ 100 µl Binding Buffer, collecting each 1 ml.

Elute 10x w/ 100 µl Election Buffer, collecting each 100 µl.

Wash column 10x w/ 100 µl Election Buffer, collecting 1 ml.

Wash column 10x w/ 100 µl Binding Buffer w/o neomycin, collecting 1 ml.

Assay 20 µl / lpc in 5 ml Scintiverse.

* columns only slightly hot afterwards

Samples:	1	Control	20 µl / 20 µl 1/4 sample	17	} Unlocking 28%
	2	Column reWash	20 µl / 100 µl	18	
	3-4	Washes	20 µl / 100 µl	19-20	
	5-14	Election	20 µl / 100 µl	21-30	
	15-16	Column Washes	20 µl / 100 µl	31-32	
			Control 9%		

Repeat exactly as above except neither election buffer contains neomycin
Control - 25% Unlocking 59%

Repeat as above w/ election buffer w/o neomycin
Control 9% Unlocking 25%

Protocol #: 3 Name: 32PMNH 14:20
Region A: LL-UL= 5.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10
Region B: LL-UL=50.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10
Region C: LL-UL= 0.0- 0.0 Lcr= 0 Bkg= 0.00 %2 Sigma=0.00
e = 2.00 QIP = SIS

S#	TIME	CPMA	CPMB	B:2S%	SIS FLAG
1	2.00	11164.0	8355.50	1.55	544.96
2	2.00	21.50	10.50	43.64	451.41
3	2.00	1360.50	1064.50	4.33	641.54
4	2.00	199.00	129.50	12.43	398.68
5	2.00	372.00	258.00	8.80	529.43
6	2.00	489.50	372.50	7.33	557.94
7	2.00	349.00	242.00	9.09	510.14
8	2.00	260.50	178.50	10.59	518.90
9	2.00	120.00	82.00	15.62	432.42
10	2.00	72.50	53.50	19.33	641.54
11	2.00	53.50	30.50	25.61	427.76
12	2.00	37.00	22.50	29.81	473.56
13	2.00	39.00	19.00	32.44	273.50
14	2.00	36.50	16.00	35.36	436.39
15	2.00	26.50	12.50	40.00	334.99
16	2.00	36.00	22.50	29.81	415.56
17	2.00	11832.0	8961.00	1.49	599.53
18	2.00	20.50	8.50	48.51	124.50
19	2.00	1274.50	951.50	4.58	523.71
20	2.00	249.00	163.50	11.06	480.42
21	2.00	555.00	406.00	7.02	543.77
22	2.00	1571.50	1116.50	4.23	469.11
23	2.00	1590.00	1200.00	4.08	566.86
24	2.00	1090.00	774.50	5.08	503.81
25	2.00	462.00	326.00	7.83	512.29
26	2.00	152.00	95.50	14.47	386.64
27	2.00	131.50	95.50	14.47	521.19
28	2.00	64.00	35.50	23.74	397.77
29	2.00	55.50	38.50	22.79	535.36
30	2.00	55.00	30.00	25.82	446.85
31	2.00	37.00	19.00	32.44	497.00
32	2.00	29.00	15.50	35.92	565.03

3A
v000(2)

Protocol #: 3 Name: 32PMNH
 Region A: LL-UL= 5.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10
 Region B: LL-UL=50.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10
 Region C: LL-UL= 0.0- 0.0 Lcr= 0 Bkg= 0.00 %2 Sigma=0.00
 Time = 2.00 QIP = SIS

S#	TIME	CPMA	CPMB	B:25%	SIS FLAG
1	2.00	5477.00	4327.00	2.15	630.49
2	2.00	1087.00	714.50	37.14	346.27
3	2.00	829.00	624.00	5.66	553.51
4	2.00	602.00	470.50	6.52	631.53
5	2.00	294.50	221.00	9.51	539.05
6	2.00	152.50	112.00	13.36	572.09
7	2.00	86.00	60.50	18.18	452.53
8	2.00	51.50	31.50	25.20	464.95
9	2.00	40.00	22.00	30.15	315.86
10	2.00	30.50	15.00	36.51	318.22
11	2.00	30.00	13.50	38.49	341.19
12	2.00	24.00	14.00	37.80	378.88
13	2.00	42.00	24.50	28.57	312.21
14	2.00	4437.50	3492.50	2.39	601.46
15	2.00	36.50	24.00	28.87	460.09
16	2.00	859.50	607.50	5.74	513.46
17	2.00	197.50	150.50	11.53	569.37
18	2.00	379.50	271.50	8.58	507.38
19	2.00	1059.00	782.50	5.06	534.56
20	2.00	1211.00	908.00	4.69	557.88
21	2.00	1002.50	705.00	5.33	524.97
22	2.00	446.50	330.50	7.78	501.70
23	2.00	254.50	180.50	10.53	542.42
24	2.00	166.50	120.50	12.88	492.37
25	2.00	113.50	79.00	15.91	501.40
26	2.00	62.50	40.00	22.36	448.20
27	2.00	53.00	35.00	23.90	487.85
28	2.00	34.00	18.50	32.88	398.72
29	2.00	66.50	41.50	21.95	491.40

38
002(3)

Protocol #: 3 Name: 32PMNH

09:44

Region A: LL-UL= 5.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10

Region B: LL-UL=50.0-1700 Lcr= 0 Bkg= 0.00 %2 Sigma=0.10

Region C: LL-UL= 0.0- 0.0 Lcr= 0 Bkg= 0.00 %2 Sigma=0.00

Time = 2.00 QIP = SIS

S#	TIME	CPMA	CPMB	B:2S%	SIS	FLAG
1	2.00	8881.00	7140.00	1.67	685.52	
2	2.00	25.00	10.00	44.72	424.10	
3	2.00	1114.00	818.50	4.94	542.31	
4	2.00	218.50	163.00	11.08	529.52	
5	2.00	303.00	243.50	8.06	625.26	
6	2.00	596.00	472.50	8.51	675.09	
7	2.00	311.00	227.00	9.38	551.88	
8	2.00	142.50	95.00	14.51	443.71	
9	2.00	91.00	62.00	17.96	440.44	
10	2.00	40.50	23.50	29.17	345.52	
11	2.00	41.00	21.50	30.50	332.35	
12	2.00	32.00	19.50	32.03	360.78	
13	2.00	25.50	12.00	40.82	552.16	
14	2.00	42.00	22.50	29.81	351.20	
15	2.00	28.50	13.00	39.22	532.13	
16	2.00	39.50	20.50	31.23	425.34	
17	2.00	1554.00	1151.00	4.17	559.42	
18	2.00	47.00	31.50	25.20	165.14	
19	2.00	114.50	85.50	15.29	537.42	
20	2.00	26.50	16.50	34.82	535.61	
21	2.00	52.50	31.00	25.40	529.57	
22	2.00	66.50	47.00	20.63	514.24	
23	2.00	98.00	71.00	16.78	565.60	
24	2.00	126.00	83.00	15.52	490.09	
25	2.00	130.00	90.50	14.87	572.18	
26	2.00	132.50	84.00	15.43	448.77	
27	2.00	132.00	83.50	15.48	407.11	
28	2.00	126.50	75.00	16.33	420.04	
29	2.00	115.50	79.50	15.86	521.63	
30	2.00	118.00	80.00	15.81	553.93	
31	2.00	76.00	53.50	19.33	597.11	
32	2.00	68.00	44.00	21.32	510.19	

Serial No.: 10/809,886

Docket No.: 19000.0058/P058

EXHIBIT D

TITLE: Native PAGE - Neomycin Effects on CLAMP Binding to Affinity Gels	DATE: 09/28/06
PURPOSE: To test the effects of different neomycin concentrations on CLAMP 609 binding to affinity gels	EXP.# BK041

REAGENTS

³²P-Neo-uu-ATP-CLAMP 609 (~ 2pmol/ul)/ BK040
 (gcuuaauacgacucacuaauaggccugggcgagaaguuuaggcc uu
 gguugggaagaaacuguggcacuucggugccagcaaccc)
 Neomycin Sepharose from BK001
 ATP Agarose from Sigma No. A2767
 Neomycin from Sigma No. N-1876 Lot 90K0854 MW 908.9
 Binding Buffer: 50mM Tris, 250mM NaCl, 5mM MgCl₂, pH 7.6
 Neomycin Solution: 10mM neomycin in Binding Buffer
 Nondenaturing gel 5x loading dye (bromophenol blue and xylene cyanol FF dyes)
 Nondenaturing 12% acrylamide gel, TBE Buffer, PAGE apparatus, Saran Wrap
 Whatman Chromatography Paper No. 3003917, Gel Drying Film from Promega No. V7131
 BioRad Model 583 Gel Dryer, Molecular Dynamics Storage Phosphor Screen & Cassette
 Molecular Dynamics PhosphorImager Model 400A (2127 MBB)

STEP	SOLUTION			VOL (ml)	TEMP (°C)	Time Period
Prepare Samples	~ 4pmol ³² P-Neo/ATP-CLAMP 50mM Tris 250mM NaCl 5mM MgCl ₂ pH 7.6	± neomycin (0.1, 0.5, 1, 10 mM)	± neomycin-sepharose ----- ± ATP-agarose	8 ul	23	10min
Load Gel	Add loading dye to above and load samples on gel, skipping a well after sample 1 to mark lane 1.			10ul	23	5min
Run Gel	Run native PAGE electrophoresis at 30mA/gel			1 gel	23	1.5hr
Dry Gel	Dry gel between layers of gel drying film on Whatman paper.			1 gel	80	1hr
Expose to Screen	Wrap dried gel in plastic wrap and place in phosphorimager cassette overnight.			1 gel	23	18hr
Image	Record image on phosphorimager. (Gel# 17010)					

DISCUSSION:

The image of this gel had background spots and the bands were less even, but affinity gel bound CLAMP was still visible retained in the wells. More CLAMP was visible in the free bands in the gel, possibly due to the lower ligand concentration of the affinity gels as compared to the derivatized microspheres. The most interesting results from this gel is seen in samples 7-11, where as neomycin concentration is increased, the bands of free CLAMP decrease in intensity, suggesting that ATP agarose binding increases. This reinforces the results from column affinity studies (BK036 & BK038) that showed the presence of free neomycin increasing the adenosine binding capacity of CLAMP 609.

C:\DATA\BK041B.GEL , Range = 0.01-900.00 Counts, 1.00x

Native PAGE - CLAMP609 binding to affinity gels
(BK041B)

12% acrylamide gel

30% acrylamide

DEPC water

10x TBE

10% ammonium persulfate

TEMED

8.0 ml

9.86 ml

2.0 ml

0.14 ml

7 μ lBinding buffer: 50mM Tris, 250mM NaCl, 5mM MgCl₂

pH 7.6

Torr's Nondenaturing Gel 5x Loading Dye

w/ bromophenol blue & xylene cyanol FF

70bp

20bp

#17010


Sample#	CLAMP609	5x dye	Neo. Seph.	Alen. Agar.	B.B.	B.B. w/ 50mM Neo.	B.B. w/ 2.5mM Neo.	final [neonup]
1	2 μ l	2 μ l	—	—	6	—	—	0mM
2	2	2	2	—	4	—	—	0
3	2	2	2	—	3.6	0.4	0.4	0.1
4	2	2	2	—	3.6	0.4	2	0.5
5	2	2	2	—	3.8	0.2	—	1
6	2	2	2	—	2	2	—	10
7	2	2	—	2	4	—	—	0
8	2	2	—	2	3.6	0.4	0.4	0.1
9	2	2	—	2	3.6	0.4	2	0.5
10	2	2	—	2	3.8	0.2	—	1
11	2	2	—	2	2	2	—	10
12	2	2	2	2	2	—	—	0
13	2	2	2	2	1.6	0.4	0.4	0.1
14	2	2	2	2	1.6	0.4	2	0.5
15	2	2	2	2	1.8	0.2	—	1
16	2	2	2	2	—	2	—	10

Allow gel to polymerize 1hr, run gel @ 30mA for 1.5 hr. Dry w/ plastic membrane 1hr @ 80°C.
Expose to phosphorimager screen overnight (18hr).

⇒ overloaded, spotty gel

\\DATA\BK041B.GEL , Range = 0.01-900.00 Counts, 1.00x

Gel #17010



**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☐ BLACK BORDERS

☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES

☐ FADED TEXT OR DRAWING

☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING

☐ SKEWED/SLANTED IMAGES

☒ COLOR OR BLACK AND WHITE PHOTOGRAPHS

☐ GRAY SCALE DOCUMENTS

☒ LINES OR MARKS ON ORIGINAL DOCUMENT

☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY

☒ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.